

PRODUCT INFORMATION

Common Name	IMGN-529, IMGN529, K7153A-SMCC-DM1
Synonyms	CD37;TSPAN26;Tspan-26
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	N/A
Host Species	Chimeric
IgG type	IgG1
Reactivity	Human
Target	CD37
Uniprot ID	P11049
Description	PE-conjugated Anti-CD37 (naratuximab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

