

PRODUCT INFORMATION

Common Name	hu38SB19;SAR-650984;SAR650984;Sarclisa
Synonyms	CD38;ADPRC1;cADPR1;ADPRC 1;T10;cADPr hydrolase 1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CD38
Uniprot ID	P28907
Description	PE-conjugated Anti-CD38(isatuximab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

