

PRODUCT INFORMATION

Common Name	TAK-079;TAK079;TAK 079;TAK-79;TAK79;TAK 79
Synonyms	CD38;ADPRC1;cADPR1;ADPRC 1;T10;cADPr hydrolase 1
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CD38
Uniprot ID	P28907
Description	PE-conjugated Anti-CD38(mezagitamab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

