

PRODUCT INFORMATION

Common Name	SPV-T3a;SPVT3a
Synonyms	T3E;TCRE;IMD18;CD3epsilon
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG2b
Reactivity	Human
Target	CD3E
Uniprot ID	P07766
Description	PE-conjugated Anti-CD3E(dafsolimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

