

PRODUCT INFORMATION

Common Name	NI-0401;NI0401;NI-401;NI401;TZLS-401;TZLS401
Synonyms	T3E;TCRE;IMD18;CD3epsilon
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CD3E
Uniprot ID	P07766
Description	PE-conjugated Anti-CD3E(foralumab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

