

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | BG-9588, hu5c8 |
| Synonyms | CD40LG;CD154;CD40L;HIGM1;IGM;IMD3;T-BAM;TNFSF5;TRAP;gp39 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Humanized |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | CD40 Ligand |
| Uniprot ID | P29965 |
| Description | PE-conjugated Anti-CD40 Ligand (ruplizumab biosimilar) mAb |
| Delivery | 3-4 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

DIMABIO CONFIDENTIAL

