

PRODUCT INFORMATION

Common Name	KPL404;KPL-404;KPL 404
Synonyms	CD40;p50;Bp50;CDW40;TNFRSF5
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	CD40
Uniprot ID	P25942
Description	PE-conjugated Anti-CD40(abiprubart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

