

PRODUCT INFORMATION

Common Name	APX 005;APX005;APX 005;APX5;APX-5;APX 5;EPI-0050;EPI0050;EPI 0050;EPI 50;EPI-50;EPI50
Synonyms	CD40;p50;Bp50;CDW40;TNFRSF5
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CD40
Uniprot ID	P25942
Description	PE-conjugated Anti-CD40(sotigalimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

