

PRODUCT INFORMATION

Common Name	Iomab-B
Synonyms	B220;CD45;CD45R;GP180;L-CA;LCA;LY5;T200
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Mus musculus
IgG type	Mouse IgG1
Reactivity	Human
Target	CD45
Uniprot ID	P08575
Description	PE-conjugated Anti-CD45 (apamistamab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

