

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | Iomab-B |
| Synonyms | B220;CD45;CD45R;GP180;L-CA;LCA;LY5;T200 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | N/A |
| Host Species | Mus musculus |
| IgG type | Mouse IgG1 |
| Reactivity | Human |
| Target | CD45 |
| Uniprot ID | P08575 |
| Description | PE-conjugated Anti-CD45 (apamistamab biosimilar) mAb |
| Delivery | 3-4 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

DIMABIO CONFIDENTIAL

