

PRODUCT INFORMATION

Common Name	Hu5F9-G4
Synonyms	CD47;MER6;IAP
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	IgG1
Reactivity	Human
Target	CD47
Uniprot ID	Q08722
Description	PE-conjugated Anti-CD47 (magrolimab biosimilar;IgG1) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

