

PRODUCT INFORMATION

Common Name	AK-117;AK117;AK 117
Synonyms	CD47;IAP;OA3;MER6
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	CD47
Uniprot ID	Q08722
Description	PE-conjugated Anti-CD47(ligufalimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

