

PRODUCT INFORMATION

Common Name	UB-421;UB421;UB 421
Synonyms	T4;IMD79;Leu-3;OKT4D;CD4mut
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CD4
Uniprot ID	P01730
Description	PE-conjugated Anti-CD4(semzuvolimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

