

PRODUCT INFORMATION

| | |
|---|--|
| Common Name | ARGX-110;ARGX110;JNJ 74494550;JNJ74494550;JNJ-74494550;JNJ-4550;JNJ4550 |
| Synonyms | CD27L;LPFS3;CD27-L;CD27LG;TNFSF7;TNLG8A |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | N/A |
| Host Species | Whole mAb |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | CD70 |
| Uniprot ID | P32970 |
| Description | PE-conjugated Anti-CD70(cusatuzumab biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

DIMABIO CONFIDENTIAL

