

PRODUCT INFORMATION

Clone ID	1E9
Target	CD72
Synonyms	Lyb-2
Host Species	Rabbit
Description	PE-conjugated Anti-CD72 antibody(1E9), IgG1 Chimeric mAb
Delivery	3-4 weeks
Uniprot ID	P21854
IgG type	Rabbit/Human Fc chimeric IgG1
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Plays a role in B-cell proliferation and differentiation.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	PE-conjugated

