

PRODUCT INFORMATION

Common Name	MEDI-9447, MEDI9447
Synonyms	NT5E, 5'-Nucleotidase, 5'-NT, NT5, NTE
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid □ PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	IgG1(K97R,L117F,L118E,P214S)
Reactivity	Human
Target	CD73
Uniprot ID	P21589
Description	PE-conjugated Anti-CD73(oleclumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

