

PRODUCT INFORMATION

Common Name	TJ-004309, TJC-4309, TJD-5
Synonyms	NT5E;5'-Nucleotidase;5'-NT;NT5;NTE
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric Humanized
IgG type	IgG1
Reactivity	Human
Target	CD73
Uniprot ID	P21589
Description	PE-conjugated Anti-CD73(uliledlimab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

