

PRODUCT INFORMATION

Common Name	DS-6000A,DS6000A
Synonyms	CADH6;K-cadherin
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	N/A
Host Species	Humanized
IgG type	IgG1(K97R)
Reactivity	Human
Target	CDH6
Uniprot ID	P55285
Description	PE-conjugated Anti-CDH6(raludotatug biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

