

PRODUCT INFORMATION

Common Name	IMMU-130, hMN-14-SN-38, hMN-14-SN-38 ADC, hMN14-CL-SN-38, hMN14-SN-38
Synonyms	CEACAM-5;CD66e;CEA;Meconium antigen 100
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	CEACAM5
Uniprot ID	P06731
Description	PE-conjugated Anti-CEACAM5 (labetuzumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

