

PRODUCT INFORMATION

Common Name	M9140;M-9140;M 9140
Synonyms	CEACAM-5;CEA;CD66e;Meconium antigen 100
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CEACAM5
Uniprot ID	P06731
Description	PE-conjugated Anti-CEACAM5(precentabart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

