

PRODUCT INFORMATION

Common Name	BIIB 059, BIIB-059, BIIB059
Synonyms	BDCA-2;CD303
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	N/A
Host Species	Humanized
IgG type	IgG1(E239D, M241L)
Reactivity	Human
Target	CLEC4C
Uniprot ID	Q8WTT0
Description	PE-conjugated Anti-CLEC4C (litifilimab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

