

PRODUCT INFORMATION

Common Name	BMS-901608, HuLuc63, PDL-063, PDL063
Synonyms	CS1;CD319;SLAMF7;CRACC
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	CS1
Uniprot ID	Q9NQ25
Description	PE-conjugated Anti-CS1 (elotuzumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

