

PRODUCT INFORMATION

Common Name	BNT-316;BNT316;BNT 316;ONC-392;ONC392;ONC 392
Synonyms	CD;GSE;GRD4;ALP55;CD152;CTLA-4;IDDM12;CELIAC3
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CTLA-4
Uniprot ID	P16410
Description	PE-conjugated Anti-CTLA-4(gotistobart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

