

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | PSB105;PSB-105;PSB-205 (cocktail with Tuvonralimab);PSB205 (cocktail with Tuvonralimab);QL1706 (cocktail with Tuvonralimab) |
| Synonyms | CD;GSE;GRD4;ALPS5;CD152;CTLA-4;IDDM12;CELIAC3 |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | N/A |
| Host Species | Whole mAb |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | CTLA-4 |
| Uniprot ID | P16410 |
| Description | PE-conjugated Anti-CTLA-4(tuvonralimab biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

DIMABIO CONFIDENTIAL

