

**PRODUCT INFORMATION**

<b>Common Name</b>	PSB105;PSB-105;PSB-205 (cocktail with Tuvonralimab);PSB205 (cocktail with Tuvonralimab);QL1706 (cocktail with Tuvonralimab)
<b>Synonyms</b>	CD;GSE;GRD4;ALPS5;CD152;CTLA-4;IDDM12;CELIAC3
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Liquid□PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	Whole mAb
<b>IgG type</b>	IgG1
<b>Reactivity</b>	Human
<b>Target</b>	CTLA-4
<b>Uniprot ID</b>	P16410
<b>Description</b>	PE-conjugated Anti-CTLA-4(tuvonralimab biosimilar) mAb
<b>Delivery</b>	6-8 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

DIMABIO CONFIDENTIAL

