

PRODUCT INFORMATION

Common Name	BMS-986253;BMS986253;HuMax-IL8;HuMax-Inflam;MDX018;MDX-018
Synonyms	IL8;NAF;GCP1;LECT;LUCT;NAP1;GCP-1;LYNAP;MDNCF;MONAP;NAP-1;SCYB8
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CXCL8
Uniprot ID	P10145
Description	PE-conjugated Anti-CXCL8(adakitug biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

