

PRODUCT INFORMATION

| | |
|-----------------------------------------|--------------------------------------------------------------------------------------------------------------|
| Common Name | DS-1024;DS1024;S-95026;S95206;Sym-004 (Combination with Futuximab);Sym004 (Combination with Futuximab) |
| Synonyms | ERBB;ERRP;HER1;mENA;ERBB1;NNCIS;PIG61;NISBD2 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Whole mAb |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | EGFR |
| Uniprot ID | P00533 |
| Description | PE-conjugated Anti-EGFR(modotuximab biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

DIMABIO CONFIDENTIAL

