

PRODUCT INFORMATION

Common Name	VB4-845
Synonyms	GA733-2, M1S2, M4S1, MIC18, TACSTD1, TROP1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	N/A
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	EPCAM
Uniprot ID	P16422
Description	PE-conjugated Anti-EPCAM(oportuzumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

