

PRODUCT INFORMATION

Common Name	KB-004, KB004
Synonyms	TYRO4;HEK4;ETK1;ETK;EK4;HEK
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	EPHA3
Uniprot ID	P29320
Description	PE-conjugated Anti-EPHA3 (ifabotuzumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

