

PRODUCT INFORMATION

Common Name	LY-3016859;LY 3016859;LY3016859
Synonyms	ER;Ep;EPR
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	EREG
Uniprot ID	O14944
Description	PE-conjugated Anti-EREG(fepixnebart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

