

PRODUCT INFORMATION

Common Name	AB-023;AB 023;AB023;AB022;AB 022;AB-022;Xisomab 3G3
Synonyms	FXI;PTA
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	F11
Uniprot ID	P03951
Description	PE-conjugated Anti-F11(gruticibart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

