

PRODUCT INFORMATION

Clone ID	19H8
Target	FGF21
Synonyms	Fibroblast growth factor 21;FGF-21
Host Species	Rabbit
Description	PE-conjugated Anti-FGF21 antibody(19H8), IgG1 Chimeric mAb
Delivery	3-4 weeks
Uniprot ID	Q9NSA1
IgG type	Rabbit/Human Fc chimeric IgG1
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	This gene encodes a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities and are involved in a variety of biological processes. This protein is a secreted endocrine factor that functions as a major metabolic regulator. The encoded protein stimulates the uptake of glucose in adipose tissue. [provided by RefSeq, Mar 2016]
Usage	Research use only
Conjugate	PE-conjugated

