

PRODUCT INFORMATION

Common Name	AMG-477;AMG477;AMB 477;REMD-477;REMD477;REMD 477
Synonyms	GR;GCR;GRL;GCCR;GCRST
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG2
Reactivity	Human
Target	GCCR
Uniprot ID	P04150
Description	PE-conjugated Anti-GCCR(volagidemab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

