

PRODUCT INFORMATION

Common Name	REGN-1033;REGN1033;REGN 1033;SAR-391786;SAR 391786;SAR391786
Synonyms	GDF8;MSLHP;MSTN
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	GDF8
Uniprot ID	O14793
Description	PE-conjugated Anti-GDF8(trevogrumab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

