

PRODUCT INFORMATION

| | |
|---|--|
| Clone ID | 21C2 |
| Target | GFAP |
| Synonyms | ALXDRD |
| Host Species | Rabbit |
| Description | PE-conjugated Anti-GFAP(68-377) antibody(21C2), Rabbit mAb |
| Delivery | 3-4 weeks |
| Uniprot ID | P14136 |
| IgG type | Rabbit IgG |
| Clonality | Monoclonal |
| Reactivity | Human |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Purification | Purified from cell culture supernatant by affinity chromatography |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | This gene encodes one of the major intermediate filament proteins of mature astrocytes. It is used as a marker to distinguish astrocytes from other glial cells during development. Mutations in this gene cause Alexander disease, a rare disorder of astrocytes in the central nervous system. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Oct 2008] |
| Usage | Research use only |
| Conjugate | PE-conjugated |

