

PRODUCT INFORMATION

Common Name	AMG133
Synonyms	GIP-R
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	IgG1(E155C,R175C,N180G,V185C)
Reactivity	Human
Target	GIPR
Uniprot ID	P48546
Description	PE-conjugated Anti-GIPR(maridebart biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

