

PRODUCT INFORMATION

Clone ID	DM61
Target	GPRC5D
Synonyms	GPRC5D
Host Species	Rabbit
Description	PE-conjugated Anti-GPRC5D antibody(DM61); Rabbit mAb
Delivery	3-4 weeks
Uniprot ID	Q9NZD1
IgG type	Rabbit IgG
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	The protein encoded by this gene is a member of the G protein-coupled receptor family; however; the specific function of this gene has not yet been determined.
Usage	Research use only
Conjugate	PE-conjugated

