

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | 4D5-8,4D5V8,Herceptin,rhuMabHER2 |
| Synonyms | ERBB2;CD340;HER-2/neu;HER2;MLN19;NEU;NGL;TKR1 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | N/A |
| Host Species | Humanized |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | Her2 |
| Uniprot ID | P04626 |
| Description | PE-conjugated Anti-Her2 (trastuzumab biosimilar) mAb |
| Delivery | 3-4 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

DIMABIO CONFIDENTIAL

