

PRODUCT INFORMATION

Common Name	FHTR2163, RG-6147, RG6417, RO-7171009, RO7171009
Synonyms	L56;Serine protease 11
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	HTRA1
Uniprot ID	Q92743
Description	PE-conjugated Anti-HTRA1(galegenimab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

