

## PRODUCT INFORMATION

<b>Clone ID</b>	DMC367
<b>Target</b>	ICAM-1
<b>Synonyms</b>	ICAM1;BB2;CD54;P3.58
<b>Host Species</b>	Rabbit
<b>Description</b>	PE-conjugated Anti-ICAM-1 antibody(DMC367); IgG1 Chimeric mAb
<b>Delivery</b>	3-4 weeks
<b>Uniprot ID</b>	P05362
<b>IgG type</b>	Rabbit/Human Fc chimeric IgG1
<b>Clonality</b>	Monoclonal
<b>Reactivity</b>	Human
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Purification</b>	Purified from cell culture supernatant by affinity chromatography
<b>Formulation &amp; Reconstitution</b>	Liquid□PBS with 0.05% Proclin300, 1% BSA
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	This gene encodes a cell surface glycoprotein which is typically expressed on endothelial cells and cells of the immune system. It binds to integrins of type CD11a : CD18; or CD11b : CD18 and is also exploited by Rhinovirus as a receptor.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

