

## PRODUCT INFORMATION

|   |   |
|---|---|
| <b>Common Name</b>                      | ASLAN004, ASLAN-004, CSL-334  |
| <b>Synonyms</b>                         | IL13R; IL13RA; CD213a1; CT19  |
| <b>Applications</b>                     | Flow Cyt  |
| <b>Recommended Dilutions</b>            | Flow Cyt 1:100  |
| <b>Formulation &amp; Reconstitution</b> | N/A   |
| <b>Host Species</b>                     | Homo sapiens  |
| <b>IgG type</b>                         | IgG4  |
| <b>Reactivity</b>                       | Human   |
| <b>Target</b>                           | IL13RA1   |
| <b>Uniprot ID</b>                       | P78552  |
| <b>Description</b>                      | PE-conjugated Anti-IL13RA1(eblasakimab biosimilar) mAb  |
| <b>Delivery</b>                         | 3-4 weeks   |
| <b>Storage &amp; Shipping</b>           | Store at 2°C-8°C for 6 months   |
| <b>Background</b>                       | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| <b>Usage</b>                            | Research use only   |
| <b>Conjugate</b>                        | PE-conjugated   |

DIMABIO CONFIDENTIAL

