

PRODUCT INFORMATION

Common Name	HuMax-IL15, AMG 714, AMG-714, HuMax-IL-15, MAb146B7, PRV-015
Synonyms	IL-15
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	IgG1(K97R, E239D, M241L)
Reactivity	Human
Target	IL15
Uniprot ID	P40933
Description	PE-conjugated Anti-IL15(ordesekimab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

