

PRODUCT INFORMATION

Common Name	AIN457
Synonyms	CTLA-8;CTLA8;IL-17;IL-17A;IL17
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	N/A
Host Species	Homo sapiens
IgG type	IgG1
Reactivity	Human
Target	IL17A
Uniprot ID	Q16552
Description	PE-conjugated Anti-IL17A(secukinumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

