

PRODUCT INFORMATION

| | |
|-----------------------------------------|-------------------------------------------------------------------------------------------------------|
| Common Name | TBC |
| Synonyms | IL17;CTLA8;IL-17;ILA17;CTLA-8;IL-17A |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | N/A |
| Host Species | Whole mAb |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | IL17A |
| Uniprot ID | Q16552 |
| Description | PE-conjugated Anti-IL17A(turenkibart biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

