

PRODUCT INFORMATION

Common Name	GR-1501;GR1501;GR 1501
Synonyms	IL17;CTLA8;IL-17;ILA17;CTLA-8;IL-17A
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	IL17A
Uniprot ID	Q16552
Description	PE-conjugated Anti-IL17A(xeligekimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

