

PRODUCT INFORMATION

Clone ID	4G2
Target	IL1A
Synonyms	IL-1 alpha;IL-1A;IL1;IL1-ALPHA;IL1F1
Host Species	Rabbit
Description	PE-conjugated Anti-IL1A antibody(4G2); IgG1 Chimeric mAb
Delivery	3-4 weeks
Uniprot ID	P01583
IgG type	Rabbit/Human Fc chimeric IgG1
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	<p>The protein encoded by this gene is a member of the interleukin 1 cytokine family. This cytokine is a pleiotropic cytokine involved in various immune responses, inflammatory processes, and hematopoiesis. This cytokine is produced by monocytes and macrophages as a proprotein, which is proteolytically processed and released in response to cell injury, and thus induces apoptosis. This gene and eight other interleukin 1 family genes form a cytokine gene cluster on chromosome 2. It has been suggested that the polymorphism of these genes is associated with rheumatoid arthritis and Alzheimer's disease. [provided by RefSeq, Jul 2008]</p>
Usage	Research use only
Conjugate	PE-conjugated

