

PRODUCT INFORMATION

Common Name	XB20-12;XB2012;XB20 12
Synonyms	IL1;IL-1A;IL1F1;IL1-ALPHA;IL-1 alpha
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	IL1A
Uniprot ID	P01583
Description	PE-conjugated Anti-IL1A(lunlekitug biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

