

PRODUCT INFORMATION

Common Name	CAN-04, CAN04
Synonyms	IL1AP; C3orf13; IL1R3; IL-1R-3; IL-1R3
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	IL1RAP
Uniprot ID	Q9NPH3
Description	PE-conjugated Anti-IL1RAP(nadunolimab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

