

## PRODUCT INFORMATION

<b>Common Name</b>	BI-655130, BI655130, BI 655130, spesolimab-sbzo
<b>Synonyms</b>	ILRL2; IL-36R; IL1RRP2; IL-36R
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Liquid □ PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	Chimeric Humanized
<b>IgG type</b>	IgG1(K97R, L117A, L118A)
<b>Reactivity</b>	Human
<b>Target</b>	IL1RL2
<b>Uniprot ID</b>	Q9HB29
<b>Description</b>	PE-conjugated Anti-IL1RL2(spesolimab biosimilar) mAb
<b>Delivery</b>	3-4 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

DIMABIO CONFIDENTIAL

