

PRODUCT INFORMATION

Common Name	Ilumetri;ILUMYA;MK-3222;MK3222;MK3222;SCH-900222;SCH900222;SCH900222;SUNPG1622;SUNPG1622;SUNPG1622;SUNPG1623;SUNPG-1622;SUNPG1623;SUNPG-1623;SUNPG1623
Synonyms	P19;SGRF;IL-23;IL-23A;IL23P19
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	IL23A
Uniprot ID	Q9NPF7
Description	PE-conjugated Anti-IL23A(tildrakizumab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

